

SOFT SHIELD® Collagen Shield by OASIS® Medical



**Corneal Comfort & Protection**

NOW CE APPROVED



- Ocular surface protection following surgery, injection, traumatic and non-traumatic corneal conditions
- Comfort of a thin, clear, and pliable collagen barrier
- Hydrates to 70% water content
- Protects the surface of the eye
- Lubricates the ocular surface with collagen while it gradually dissolves

Available in Quick Shield, 12 hour, 24 hour, and 72 hour dissolution times

Diameter: 14.5mm    Compound base curve: 9mm

OASIS  True Surgical Solutions.