

3-MONTHLY REPLACEMENT

SAPHIR

SILICONE HYDROGEL



Saphir provides comfortable, healthy contact lens wear to patients accustomed to a 3-monthly lens replacement.


-  **SPHERIC**
-  **MULTIFOCAL**
-  **TORIC**
-  **MULTIFOCAL TORIC**

PARAMETERS

Base curves (mm)	6.80 to 9.80 (0.30)
Diameters (mm)	13.00 to 16.00 (0.50)
Spheres (D)	±30.00 (0.25)
Cylinders	-0.75 to -8.00 (0.25)
Axes (°)	All (5°)
Additions	0.50 to 4.00 (0.50) CD/CN

MATERIAL





Type	Filcon 5B (60) [75%]
DK (ISO 9913-1-1998)	60
DK/t (-3.00 D)	50
Water Content	75%
Central Thickness (-3.00 D)	0.12
Cof	0.04
Modulus	0.29
Handling tint	No
Pack size	Single and 2-pack
Manufacturing Process	Lathed



Need fitting advice?

Check the Fitting Guide section of the catalogue.

OPTICAL DESIGN



Lens design parameters may change depending on the power

CONTACT YOUR BUSINESS DEVELOPMENT MANAGER FOR DETAILS

Customer Care: T. 0800 328 0610 F. 0800 328 0649
International Customer Care: T. +34 916 496 123
mkservices@markennovy.com

markennovy