

# saline



- 1 *It is a sterile, isotonic and buffered pH solution*
- 2 *To rinse, hydrate and preserve all types of contact lenses. It is also used for weekly enzymatic cleaning with Avizor enzyme.*
- 3 *It is a very useful supplementary solution for users of all types of contact lenses as it eliminates, for example, any dust particles.*

## Technical Data:

Avizor Saline Solution contains a combination of elements that make it effective:

- **Cloruro Sódico:**  
for tonicity and osmolarity.
- **A tampon solution:**  
which controls the solution's pH so that it adapts to that of the eye.
- **Polyhexanide:**  
which is a preservative (0,00005%) (only in multidose format).
- **Unidose format is PRESERVATIVE FREE.**



## Formatos

Box of 24 Units - 350 ml  
Box of 30 Units - 500 ml  
Box of 16 Units - 30x5 ml