

4 F

A popular high-quality custom made RGP lens, based on a 4-curve design with an aspheric periphery.

CAN AID THE FOLLOWING

- Myopia
- Hyperopia
- Astigmatism

EXTRAS AVAILABLE

- Toric Flange
- Back surface Toric
- Front surface Toric
- Bi-Toric

DESIGN

- Spherical blended front surface
- Aspheric back surface

MATERIAL OPTIONS

- Optimum Classic Roflufocon A 3 (26)
- Optimum Comfort Roflufocon C 3 (65)
 - Optimum Extra Roflufocon D 3 (100)
- Optimum Extreme Roflufocon E 3 (125)
 - Optimum Infinite Tisilfocon A 3 (180)

PARAMETERS

- Base Curve 7.05, 7.12, 7.20, 7.27, 7.35, 7.42, 7.50, 7.57, 7.65, 7.72, 7.80, 7.87, 7.95, 8.02, 8.10, 8.17, 8.25, 8.32, 8.40, 8.47
- Diameter 8.60 to 10.60 (0.10 steps)
- Power $-/+ 30.00D$ (0.25D steps) - higher powers may be available.
- Handling tints Blue, Clear, Green, Grey (material dependent)

All measurements are in mm unless otherwise stated.

FITTING

Refer to our RGP fitting guide available on our website.

MODALITY

12 months maximum. Replacement frequency is at the discretion of the eye care professional.

